
AESFULAP[®] MXE[®]

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- (1) Fink U. Plasmapore: A plasma-sprayed microporous titanium coating to improve the long-term stability. *Actualités en Biomateriaux*. 1996;III:97-104.
- (2) Swamy G, Pace A, Quah C, Howard P. The Bicontact cementless primary total hip arthroplasty: long-term results. *Int Orthop*. 2012 May;36(5):915-20.
- (3) Simank HG, Stuber M, Frahm R, Helbig L, van Lenthe H, Müller R. The influence of surface coatings of dicalcium phosphate (DCPD) and growth and differentiation factor-5 (GDF-5) on the stability of titanium implants in vivo. *Biomaterials*. 2006 Jul;27(21):3988-94.
- (4) Pilz M, Staats K, Tobudic S, Assadian O, Prestler E, Windhager R, Holinka J. Zirconium Nitride Coating Reduced Staphylococcus epidermidis Biofilm Formation on Orthopaedic Implant Surfaces: An In Vitro Study. *Clin Orthop Relat Res*. 2019 Feb;477(2):461-466.
- (5) Reich J, Hovy L, Lindenmaier HL, Zeller R, Schwiesau J, Thomas P, Grupp TM. Präklinische Ergebnisse beschichteter Knieimplantate für Allergiker [Preclinical evaluation of coated knee implants for allergic patients]. *Orthopade*. 2010 May;39(5):495-502. German.
- (6) Grupp TM, Schroeder C, Kyun Kim T, Miehke RK, Fritz B, Jansson V, Utzschneider S. Biotribology of a mobile bearing posterior stabilised knee design-effect of motion restraint on wear, tibio-femoral kinematics and particles. *J Biomech*. 2014 Jul 18;47(10):2415-23.
- (7) Puente Reyna AL, Fritz B, Schwiesau J, Schilling C, Summer B, Thomas P, Grupp TM. Metal ion release barrier function and biotribological evaluation of a zirconium nitride multilayer coated knee implant under highly demanding activities wear simulation. *J Biomech*. 2018 Oct 5;79:88-96.
- (8) Lütznert J, Beyer F, Lütznert C, Tille E, Postler AE. A Novel Multilayer-Coating for Total Knee Arthroplasty Implants is Safe - 10-Year Results From a Randomized-Controlled Trial. *J Arthroplasty*. 2023 Jan;38(1):90-95.e1.
- (9) Thomas P, Hisgen P, Kiefer H, Schmerwitz U, Ottersbach A, Albrecht D, Summer B, Schinkel C. Blood cytokine pattern and clinical outcome in knee arthroplasty patients: comparative analysis 5 years after standard versus "hypoallergenic" surface coated prosthesis implantation. *Acta Orthop*. 2018 Dec;89(6):646-651.
- (10) Lütznert J, Beyer F, Postler A, Summer S, Lütznert C. 7 Year Results Of A Novel Coating In TKA. Poster #161 EFORT 2021.
- (11) AOANJRR annual report 2021.
- (12) Utzschneider S, Becker F, Grupp TM, Sievers B, Paulus A, Gottschalk O, Jansson V. Inflammatory response against different carbon fiber-reinforced PEEK wear particles compared with UHMWPE in vivo. *Acta Biomater*. 2010 Nov;6(11):4296-304.
- (13) Grupp TM, Giurea A, Miehke RK, Hintner M, Gaisser M, Schilling C, Schwiesau J, Kaddick C. Biotribology of a new bearing material combination in a rotating hinge knee articulation. *Acta Biomater*. 2013 Jun;9(6):7054-63.
- (14) Grupp TM, Holderied M, Mulliez MA, Streller R, Jäger M, Blömer W, Utzschneider S. Biotribology of a vitamin E-stabilized polyethylene for hip arthroplasty - Influence of artificial ageing and third-body particles on wear. *Acta Biomater*. 2014 Jul;10(7):3068-78.
- (15) Baghdadi J, Alkhateeb S, Roth A; VITAS-Group; Jäger M. Cup positioning and its effect on polyethylene wear of vitamin E- and non-vitamin E-supplemented liners in total hip arthroplasty: radiographic outcome at 5-year follow-up. *Arch Orthop Trauma Surg*. 2022 Apr 10.
- (16) Puente Reyna AL, Schwiesau J, Altermann B, Grupp TM. Mechanical properties and wear resistance of a vitamin E blended and moderately cross-linked polyethylene for total knee arthroplasty. Poster IORS 2022.
- (17) de Steiger R, Lorimer M, Graves SE. Cross-Linked Polyethylene for Total Hip Arthroplasty Markedly Reduces Revision Surgery at 16 Years. *J Bone Joint Surg Am*. 2018 Aug 1;100(15):1281-1288.
- (18) Diamond OJ, Howard L, Masri B. Five cases of tibial post fracture in posterior stabilized total knee arthroplasty using Prolong highly cross-linked polyethylene. *Knee*. 2018 Aug;25(4):657-662.
- (19) MacDonald DW, Higgs GB, Chen AF, Malkani AL, Mont MA, Kurtz SM. Oxidation, Damage Mechanisms, and Reasons for Revision of Sequentially Annealed Highly Crosslinked Polyethylene in Total Knee Arthroplasty. *J Arthroplasty*. 2018 Apr;33(4):1235-1241.
- (20) MacDonald D, Hanzlik J, Sharkey P, Parvizi J, Kurtz SM. In vivo oxidation and surface damage in retrieved ethylene oxide-sterilized total knee arthroplasties. *Clin Orthop Relat Res*. 2012 Jul;470(7):1826-33.
- (21) Kop AM, Pabbruwe MB, Keogh C, Swarts E. Oxidation of Second Generation Sequentially Irradiated and Annealed Highly Cross-Linked X3™ Polyethylene Tibial Bearings. *J Arthroplasty*. 2015 Oct;30(10):1842-6.
- (22) Rowell SL, Muratoglu OK. Investigation of surgically retrieved, vitamin E-stabilized, crosslinked UHMWPE implants after short-term in vivo service. *J Biomed Mater Res B Appl Biomater*. 2016 Aug;104(6):1132-40.
- (23) Muratoglu OK, Bragdon CR, O'Connor DO, Perinchieff RS, Jasty M, Harris WH. Aggressive wear testing of a cross-linked polyethylene in total knee arthroplasty. *Clin Orthop Relat Res*. 2002 Nov;404(4):89-95.
- (24) Orita K, Minoda Y, Sugama R, Ohta Y, Ueyama H, Takemura S, Nakamura H. Vitamin E-infused highly cross-linked polyethylene did not reduce the number of in vivo wear particles in total knee arthroplasty. *Bone Joint J*. 2020 Nov;102-B(11):1527-1534.
- (25) Matthews JB, Besong AA, Green TR, Stone MH, Wroblewski BM, Fisher J, Ingham E. Evaluation of the response of primary human peripheral blood mononuclear phagocytes to challenge with in vitro generated clinically relevant UHMWPE particles of known size and dose. *J Biomed Mater Res*. 2000 Nov;52(2):296-307.
- (26) Iwakiri K, Minoda Y, Kobayashi A, Sugama R, Iwaki H, Inori F, Hashimoto Y, Ohashi H, Ohta Y, Fukunaga K, Takaoka K. In vivo comparison of wear particles between highly crosslinked polyethylene and conventional polyethylene in the same design of total knee arthroplasties. *J Biomed Mater Res B Appl Biomater*. 2009 Nov;91(2):799-804.
- (27) Schwiesau J, Fritz B, Bergmann G, Puente Reyna AL, Schilling C, Grupp TM. Influence of radiation conditions on the wear behaviour of Vitamin E treated UHMWPE gliding components for total knee arthroplasty after extended artificial aging and simulated daily patient activities. *J Mech Behav Biomed Mater*. 2021 Oct;122:104652.
- (28) Grupp TM, Fritz B, Kutzner I, Schilling C, Bergmann G, Schwiesau J. Vitamin E stabilised polyethylene for total knee arthroplasty evaluated under highly demanding activities wear simulation. *Acta Biomater*. 2017 Jan 15;48:415-422.
- (29) Schwiesau J, Fritz B, Kutzner I, Bergmann G, Grupp TM. CR TKA UHMWPE wear tested after artificial aging of the vitamin E treated gliding component by simulating daily patient activities. *Biomed Res Int*. 2014;2014:567374.
- (30) Lerf R, Zurbrugg D, Delfosse D. Use of vitamin E to protect cross-linked UHMWPE from oxidation. *Biomaterials*. 2010 May;31(13):3643-8.
- (31) Schwiesau J, Schilling C, Kaddick C, Utzschneider S, Jansson V, Fritz B, Blömer W, Grupp TM. Definition and evaluation of testing scenarios for knee wear simulation under conditions of highly demanding daily activities. *Med Eng Phys*. 2013 May;35(5):591-600.
- (32) Jasty M, Smith E. Wear particles of total joint replacements and their role in periprosthetic osteolysis. *Curr Opin Rheumatol*. 1992 Apr;4(2):204-9.
- (33) Data on file, Aesculap Lab Test report T106.
- (34) Schwiesau J, Fritzb, Kutzner I, Bergmann G, Grupp TM. Wear mode of artificially aged UHMWPE knee bearings Evaluated by The Simulation of daily Patient activities. #0084 ESB 2015.
- (35) Data on file, Aesculap Lab Test report T116.
- (36) Zietz C, Reinders J, Schwiesau J, Paulus A, Kretzer JP, Grupp T, Utzschneider S, Bader R. Experimental testing of total knee replacements with UHMWPE inserts: impact of severe wear test conditions. *J Mater Sci Mater Med*. 2015 Mar;26(3):134.

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B. Braun Deutschland GmbH & Co. KG | Tel.: (0 56 61) 9147-70 00 | E-Mail: info.de@bbraun.com | www.bbraun.de
Aesculap AG | Am Aesculap-Platz | 78532 Tuttlingen | Deutschland