

Contact: Philipp Deist
Fon: +49 5661 71-4827
Fax: +49 5661 75-4827
Email: philipp.deist@bbraun.com
Internet: <http://www.bbraun.com>

Date: May 13, 2016

Travel-Certificate for the IryPump® bowel irrigation system

The IryPump® S Set is a medical device. Designed for ostomy patients to irrigate their bowel. The ostomy-irrigation is a therapy-concept for patients to get back their quality of life and their control of stool emptying.

For a successful therapy, the owner has to use the IryPump® S Set every or every second day. Also, when he is on travelling.

The IryPump® S Set contains:

- 1.) The electrical IryPump® S station with a fixed Li-Ion Accu (3,7 V DC, 1620mAh)
- 2.) Removable and expandable water-container
- 3.) Power supply (100-240V AC / 5,1 V DC)
- 4.) IryCone+ (cone and tubesystem)
- 5.) additional products for carrying the ostomy (e.g. stomapouches, baseplates, etc.)

Included in the IryPump® S station is an electrical water pump, which inflates the water into the bowel with a maximal water pressure of approximately 480 mbar and a maximal water flow of approximately 550 ml/min at atmospheric counter pressure, stepless adjustable.

All parts of the IryPump® S Set are CE-marked according to Directive 93/42/EEC concerning Medical Devices.

B. Braun Melsungen AG

i. V.

i. A.

Henning Ross
Head of Marketing Ambulant & Nephrologie
Marketing Germany

Philipp Deist
Marketing Manager Continence
Marketing Ambulant & Nephrologie

Chairman of Supervisory Board:
Prof. Dr. h.c. Ludwig Georg Braun

Executive Board:
Prof. Dr. Heinz-Walter Große
(Chairman)
Dr. Annette Beller
Anna Maria Braun
(Deputy Board Member)

Otto Philipp Braun
Prof. Dr. Hanns-Peter Knaebel
Dr. Meinrad Lugan
Caroll H. Neubauer, LL.M.
Markus Strotmann

Corporate Office: Melsungen
Register Court: Local Court Fritzlar
HRB 11 000
WEEE-Reg.-No. DE 42690900

Address:
B. Braun Melsungen AG
Carl-Braun-Straße 1
34212 Melsungen
Germany