



## DISPOSABLE SKIN STAPLER 35W

# DISPOSABLE SKIN STAPLER 35W



## SKIN STAPLER

- Lightweight and ergonomic
- Grip design compatible with smaller hands
- Robust assembly made with high quality materials
- Angled head to improve visibility
- Plastic cartridge that provides a direct view of remaining staples

## STAPLE REMOVERS

- Available as single use and reusable remover
- Can be used for virtually all brands of surgical skin staplers
- Intuitive easy to use device for a reliable and quick staple removal

## ORDER INFORMATION

PRODUCT DESCRIPTION	STAPLE INFORMATION TYPE/NO.	STAPLE DIMENSIONS CROWN/LEG	ARTICLE NO.	UNITS PER BOX
B. Braun Skin Stapler 35W	Wide/35	7.2 mm/4.9 mm	783106	5
Reusable Skin Stapler Remover			783102	1
Staple Remover C-W			783104	20



Manufacturer acc. to directive 93/42/EEC:  
Ningbo Advan Electrical Co., Ltd.  
Industrial Development Zone | Fuhai Town | CiXi City | China



Shanghai International Holding Corp. GmbH (Europe)  
Eiffestraße 80 | 20537 Hamburg | Germany

Distributed by:  
B. Braun Surgical, S.A. | Carretera de Terrassa, 121 | 08191 Rubí | Spain  
Phone +34 93 586 6200 | Fax +34 93 588 1096 | [www.bbraun.com](http://www.bbraun.com)

## AESCULAP<sup>®</sup> – a B. Braun brand

Aesculap AG | Am Aesculap-Platz | 78532 Tuttlingen | Germany  
Phone +49 7461 95-0 | Fax +49 7461 95-2600 | [www.aesculap.com](http://www.aesculap.com)

The main product trademark "Aesculap" is a registered trademark of Aesculap AG.

Subject to technical changes. All rights reserved. This brochure may only be used for the exclusive purpose of obtaining information about our products. Reproduction in any form partial or otherwise is not permitted.